

Disinfecting Products

Kill Claims and Contact Times

A full portfolio of actives and forms to meet all your cleaning and disinfecting needs.* Kills over 80 microorganisms, including *C. difficile* spores, Norovirus, TB and 22 antibiotic-resistant organisms!



		Bleach-Based Products				Hydrogen Peroxide-Based Products		Quat-Based Products	Alcohol-Based Products	Quat Alcohol-Based Products
		Clorox Healthcare® Bleach Germicidal Wipes	Clorox Healthcare® Bleach Germicidal Cleaners	Clorox Healthcare® Fuzion® Cleaner Disinfectant	Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectants	Clorox Healthcare® VersaSure™ Cleaner Disinfectant Wipes	Clorox Healthcare® Citrace® Hospital Disinfectant & Sanitizer	Clorox Healthcare® EZ-Kill® Quat Alcohol Cleaner Disinfectant Wipes
EPA-registered to kill:										
Bacteria	<i>Acinetobacter baumannii</i>	30 sec	1 min	1 min	1 min	1 min	30 sec		5 min	
	Multidrug-resistant <i>Acinetobacter baumannii</i>							2 min		
	<i>Bordetella bronchiseptica</i>									
	<i>Bordetella pertussis</i>	30 sec	1 min	1 min				2 min		
	<i>Burkholderia cepacia</i>			1 min				2 min		
	<i>Campylobacter jejuni</i>	30 sec	1 min	1 min		1 min	30 sec	2 min		
	<i>Chlamydia psittaci</i>									
	<i>Clostridium difficile</i> spores	3 min	3 min	2 min	3 min					
	<i>Enterobacter aerogenes</i>	30 sec	1 min	1 min	1 min	30 sec	30 sec	2 min	5 min	
	<i>Enterobacter cloacae</i>	30 sec								
	<i>Enterobacter cloacae</i> NDM-1		1 min	1 min						
	<i>Enterococcus faecalis</i>			1 min						
	Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE)	30 sec	1 min	1 min	1 min	30 sec	30 sec	2 min	5 min	2 min
	Multidrug-resistant <i>Enterococcus faecium</i>	30 sec	1 min	1 min	1 min	30 sec	30 sec	2 min	5 min	
	Vancomycin-resistant <i>Enterococcus faecium</i> (VRE)		1 min							
	<i>Enterococcus hirae</i>	30 sec	1 min	1 min	1 min					
	<i>Escherichia coli</i>	30 sec				30 sec	30 sec	2 min		
	Carbapenem-resistant <i>Escherichia coli</i>							2 min		
	<i>Escherichia coli</i> O157:H7	30 sec		1 min		30 sec	30 sec	2 min		2 min
	ESBL-producing <i>Escherichia coli</i>	30 sec	1 min		1 min	30 sec	30 sec		5 min	
	<i>Escherichia coli</i> NDM-1	30 sec	1 min			30 sec	30 sec			
	<i>Klebsiella oxytoca</i>					30 sec	30 sec			
	<i>Klebsiella pneumoniae</i>	30 sec	1 min		1 min			2 min	5 min	
	Carbapenem-resistant <i>Klebsiella pneumoniae</i>	30 sec	1 min	1 min		30 sec	30 sec	2 min		
	ESBL-producing <i>Klebsiella pneumoniae</i>	30 sec	1 min			30 sec	30 sec			
	Multidrug-resistant <i>Klebsiella pneumoniae</i>					30 sec	30 sec	2 min	5 min	
	<i>Klebsiella pneumoniae</i> NDM-1	30 sec	1 min	1 min		30 sec	30 sec			
	<i>Legionella pneumophila</i>	30 sec		1 min				2 min		
	<i>Listeria monocytogenes</i>	30 sec		1 min						
	<i>Mycobacterium bovis</i> (TB)	3 min		1 min	2 min	5 min	4 min	2 min		2 min
	<i>Proteus mirabilis</i>	30 sec	1 min	1 min						
	<i>Pseudomonas aeruginosa</i>	30 sec	1 min	1 min	1 min	30 sec	30 sec	2 min	5 min	2 min
	<i>Salmonella enterica</i> (tested as <i>Salmonella choleraesuis</i>)	30 sec	1 min	1 min	1 min	30 sec	30 sec	2 min	5 min	2 min
	<i>Serratia marcescens</i>	30 sec	1 min	1 min	1 min			2 min		
	<i>Shigella dysenteriae</i>	30 sec		1 min				2 min		
	<i>Staphylococcus aureus</i>	30 sec	1 min	1 min	1 min	1 min	1 min	2 min	5 min	2 min
	Linezolid-resistant <i>Staphylococcus aureus</i> (LRSA)	30 sec				1 min	1 min	2 min		
	Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) (ATCC 33592)	30 sec	1 min	1 min	1 min	1 min	1 min	2 min	5 min	
	Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) (ATCC 33591)					1 min	1 min			2 min
	Community-acquired Methicillin-resistant <i>Staphylococcus aureus</i> (CA-MRSA) (USA 300)	30 sec				1 min	1 min	2 min		
	Community-acquired Methicillin-resistant <i>Staphylococcus aureus</i> (CA-MRSA) (USA 400)	30 sec						2 min	5 min	
	Multidrug-resistant <i>Staphylococcus aureus</i>			1 min					5 min	
	Vancomycin-resistant <i>Staphylococcus aureus</i> (VRSA)	30 sec	1 min					2 min		
	<i>Vancomycin intermediate-resistant Staphylococcus aureus</i> (VISA)	30 sec	1 min	1 min				2 min		
	<i>Staphylococcus epidermidis</i> (Coagulase-negative staphylococci)	30 sec	1 min	1 min			1 min	2 min		
	<i>Stenotrophomonas maltophilia</i>		1 min	1 min		30 sec	30 sec			
	<i>Streptococcus equi</i> (causes distemper)									
	<i>Streptococcus pneumoniae</i>	30 sec	1 min	1 min						
	Multidrug-resistant <i>Streptococcus pneumoniae</i>			1 min						
	Penicillin-resistant <i>Streptococcus pneumoniae</i> (PRSP)					30 sec	30 sec	2 min	5 min	
	<i>Streptococcus pyogenes</i>	30 sec	1 min	1 min	1 min	30 sec	30 sec	2 min		

*Use as directed on hard, nonporous surfaces. Based on EPA Master Labels as of November 2019.

Disinfecting Products

Kill Claims and Contact Times

A full portfolio of actives and forms to meet all your cleaning and disinfecting needs.* Kills over 80 microorganisms, including *C. difficile* spores, Norovirus, TB and 22 antibiotic-resistant organisms!



		Bleach-Based Products				Hydrogen Peroxide-Based Products		Quat-Based Products	Alcohol-Based Products	Quat Alcohol-Based Products
		Clorox Healthcare® Bleach Germicidal Wipes	Clorox Healthcare® Bleach Germicidal Cleaners	Clorox Healthcare® Fuzion® Cleaner Disinfectant	Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectants	Clorox Healthcare® VersaSure™ Cleaner Disinfectant Wipes	Clorox Healthcare® Citrace® Hospital Disinfectant & Sanitizer	Clorox Healthcare® EZ-Kill® Quat Alcohol Cleaner Disinfectant Wipes
EPA-registered to kill:										
Viruses	Adenovirus Type 2	1 min	1 min	1 min	1 min					
	Avian Influenza A Virus strain H3N2									
	Avian Influenza A strain H5N1			1 min		1 min	30 sec	30 sec		
	Avian Influenza A Virus	1 min	1 min		1 min					
	Canine Distemper Virus									
	Canine Para Influenza virus (causes canine cough)									
	Canine Parvovirus	3 min	1 min	1 min	1 min					
	Coxsackie Type B3 Virus								5 min	
	Cytomegalovirus					30 sec	30 sec	30 sec		
	Enterovirus D68	1 min	1 min	1 min						
	Feline Leukemia Virus									
	Feline Panleukopenia Virus	3 min	1 min		1 min					
	Hepatitis A Virus (HAV)	1 min	1 min	1 min	1 min				5 min	
	Hepatitis B Virus (HBV) [†]	1 min	1 min	1 min	1 min	30 sec	30 sec	30 sec	1 min	2 min
	Hepatitis C Virus (HCV) [†]	1 min	1 min	1 min	1 min	30 sec	30 sec	30 sec		2 min
	Herpes Simplex Virus (HSV-1)			1 min		30 sec	30 sec	30 sec	5 min	
	Herpes Simplex Virus (HSV-2)	1 min	1 min		1 min	30 sec	30 sec	30 sec		2 min
	Human Coronavirus	1 min		1 min		30 sec	30 sec			
	Human Immunodeficiency Virus Type 1 (HIV-1)	30 sec	1 min	1 min	1 min	30 sec	30 sec	30 sec	1 min	2 min
	Infectious Bovine Rhinotracheitis Virus									
	Influenza A Virus	1 min	1 min	1 min	1 min	30 sec	30 sec	30 sec	5 min	2 min
	Influenza A Virus (H1N1)		1 min		1 min	30 sec	30 sec	30 sec		
	Influenza B Virus	1 min	1 min	1 min	1 min	30 sec	30 sec			
	Measles Virus	1 min	1 min	1 min				30 sec		
	Middle East Respiratory Syndrome Coronavirus (MERS-CoV)	1 min	1 min	1 min				30 sec		
	Mumps Virus							30 sec		
	Norovirus [‡]	1 min	1 min	1 min	1 min	3 min	1 min	5 min		
	Poliovirus Type 1 (Chat strain, ATCC VR-1562)		1 min							
	Poliovirus Type 1 (Mahoney)	1 min	1 min		1 min					
	Respiratory Syncytial Virus (RSV)	1 min		1 min		30 sec	30 sec	30 sec	5 min	
	Rhinovirus	1 min	1 min	1 min	1 min	1 min	1 min		5 min	
	Rotavirus	1 min	1 min	1 min	1 min	1 min	1 min	2 min	5 min	
	SARS-associated Coronavirus (SARS-CoV)	1 min	1 min	1 min						
Fungi	<i>Aspergillus brasiliensis</i>	5 min								
	<i>Candida albicans</i>	3 min	1 min	1 min		5 min	3 min	2 min		
	<i>Candida auris</i>					5 min	3 min			
	<i>Candida glabrata</i>	3 min		1 min						
	<i>Trichophyton mentagrophytes</i> (Athlete's Foot Fungus)	3 min	1 min	1 min	1 min	5 min	3 min		5 min	5 min
TOTALS:		60 Claims	51 Claims	53 Claims	32 Claims	43 Claims	44 Claims	45 Claims	24 Claims	13 Claims

* Use as directed on hard, nonporous surfaces.

[†] Obtained via surrogate

[‡] Murine Norovirus Strain MNV-G

Based on EPA Master Labels as of November 2019.