

Advancing
performance by design.

Duraprene™ CP

Sterile Neoprene Cleanroom Gloves

Recommended for ISO 5 /
Class 100 or higher

Features and benefits:

- Powder-free, synthetic neoprene
- Clean-processed and chlorinated
- Available in hand-specific sizes
- Proprietary interlocking, beaded cuff design reduces roll-down
- AQL of 1.0 for freedom from holes
- Sterile – gamma irradiated to SAL 10^{-6}
 - Lot specific Certificates of Conformity and Sterility are available via email at certificates@cardinalhealth.com
 - Meets ANSI/AAMI/ISO standards
- Manufactured in an ISO 13485 certified production facility
- Produced on an exclusive hand mold with an independent thumb design to allow for an anatomical fit and natural hand movement



CardinalHealth
Essential to care™



ISO-MED

ISO-MED +1 800-797-1405 | 1220 Graphite Drive, Corona, CA 92881

www.iso-med.com

Backed by the expertise of Cardinal Health, Duraprene™ CP Controlled Environment Gloves meet all relevant ASTM standards. Documentation is available upon request.

Physical dimensions and properties

Dimensions		Nominal	ASTM D 3577
Length		12 in. / 305 mm	Pass
Thickness	Finger	8.0 mils / 0.20 mm	Pass
	Palm	7.9 mils / 0.20 mm	Pass
	Cuff	6.7 mils / 0.17 mm	Pass
Properties		Nominal	ASTM D 3577
Tensile strength	Unaged	22 MPa	Pass
	Aged	18 MPa	Pass
Ultimate elongation (elasticity)	Unaged	845%	Pass
	Aged	766%	Pass
Stress at 500% elongation (modulus)	Unaged	2.8 MPa	Pass
Particle count (per cm ² > 0.5 microns)		3,251 particles/cm ²	IENT RP-CC005
Ionic extractables (mcg/cm ²)		Chemical name	IENT RP-CC005
		Bromide	0.000
		Calcium	1.364
		Chloride	0.566
		Fluoride	0.000
		Magnesium	0.002
		Nitrate	0.025
		Nitrite	0.001
		Phosphate	0.002
		Potassium	0.008
		Sodium	0.044
		Sulfate	0.191

35-month expiration from date of manufacture

Store in unopened cartons, away from excessive cold/heat, ozone and ultraviolet light.

Preferred temperature range: >40°F (4.44°C) and <100°F (37.77°C)

Catalog No.	Size
2Y72PN55	5 ½
2Y72PN60	6
2Y72PN65	6 ½
2Y72PN70	7
2Y72PN75	7 ½
2Y72PN80	8
2Y72PN85	8 ½
2Y72PN90	9

200 pairs per case, 4 polyethylene bags per case



cardinalhealth.com/cleanroomgloves

© 2017 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE and DURAPRENE are trademarks or registered trademarks of Cardinal Health. All other marks are the property of their respective owners. Lit. No. 2PER117-696366 (09/2017)

Cardinal Health
Waukegan, IL 60085