




Cotton Tipped Applicator

Wood handle



Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch handle by an aqueous based adhesive.

Dimensions

6" Overall Length (76.2 mm)

Tip

- 0.188" Dia. (4.78mm)
- 0.625" Length (15.88 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.812" Length (147.62 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Made in USA

Applications

- Applying salves or ointments
- Wound care
- Application of topical dressings
- Specimen collection
- Applying adhesives or other liquids
- Removing dirt, contamination, fluxes or resins
- Cleaning with IPA or other solvents

Ref #	Packaging
806-WC Non-sterile	10/10/100 (10 boxes per case, 100 bags per box, 100 swabs per bag)



Puritan[®]



Cotton Tipped Applicator

Wood handle

Ref #: 806-WC

NVR

NON-VOLATILE RESIDUE

Test method: IDEMA method M7-98

NVR Solvent - Hexane

NVR 12 µg/swab

The minimum detection limit for NVR is 10 µg/swab

ION CHROMATOGRAPHY TESTING RESULTS

IONIC CLEANLINESS

Test method: IDEMA M13-99 (anions) M12-99 (cations)

Anions

Fluoride	ND
Chloride	0.9000 µg/swab
Nitrate	ND
Bromide	2.5000 µg/swab
Nitrite	0.3250 µg/swab
Sulfate	4.3000 µg/swab
Phosphate	11.2500 µg/swab

Cations

Lithium	ND
Sodium	10.8000 µg/swab
Ammonium	0.7500 µg/swab
Potassium	13.7550 µg/swab
Magnesium	2.7450 µg/swab
Calcium	6.8250 µg/swab

ND - Non Detected

Test procedures available upon request.



Puritan