
























































# Clinical Application Matrix

## Exam gloves


### High Use portfolio

Product Name	Fingertip Thickness (+/-0.2mil) v	Department	Nursing Floors	Intensive Care Units (ICUs), Open Heart	OR: Surgery, Pre Op Holding, PACU, OutPatient Surgery, Anesthesia	Radiology, Nuclear	Geriatric Unit, Rehabilitation Unit, Alzheimer's Unit, Spinal Cord Injury Unit	Labs: Hematology, Microbiology, Blood Bank, Phlebotomy, Pathology	Cardiac Cath Lab, EP Labs, Interventional Radiology	Sterile Processing (Clean Side)	Obstetrics, Labor and Delivery	Emergency Departments (ED)	Housekeeping
<b>FLEXAL™ Touch</b> 88RT01XS-88RT06XXL	USP <800> 30.2 min 3.5mil												
<b>FLEXAL™ Comfort</b> 88FS01XS-88FS06XXL	USP <800> 3.7mil												
<b>FLEXAL™</b> 88TN01XS-88TN06XXL 88TNF01XS-88TNF06XXL	USP <800> 3.7mil												
<b>ESTEEM™ Comfort</b> 88ES01XS-88ES06XXL	USP <800> 50.6 min 4.5mil												
<b>ESTEEM™ Stretch</b> 88S4NXSB-88S9NXXB	USP <800> 4.5mil												
<b>ESTEEM™ Moisturize</b> N88RX01T-N88RX06T	USP <800> 32.2 min 4.7mil												
<b>ESTEEM™</b> 8895NB-8899NXXB	USP <800> 32.2 min 4.7mil												

All products are selected given the department types listed, however protection levels may vary. Glove recommendations will depend on your specific facility and protocols.

#### Key

Protection against blood and bodily fluid



USP <800> compliant



Thiotepa minimum breakthrough detection time in minutes



Fingertip texture



Fully textured



† Meets the ASTM D6978 standard























This Cardinal Health product helps healthcare facilities comply with the USP <800> standard for the safe handling of hazardous drugs.

This document was developed by a group of clinicians. It reflects current best practices of exam glove usage per hospital department.

\* Recommended glove usage for these departments are better suited for our specialty products, as they require additional and/or specific protection.

# Clinical Application Matrix – Exam gloves

## Specialty portfolio

Product Name	Fingertip Thickness (+/-0.2mil) v	Application	Sterile procedures	High risk areas (high blood and bodily fluid exposure areas)	Chemotherapy drug administration	Heavy dressing changes	Decontamination areas	Phlebotomy	Endoscopy cleaning
<b>FLEXAL™ Feel (Nitrile)</b> 88TT20XS-88TT24XL	3.3mil 								
<b>Sterile singles (Nitrile)</b> N8820-N8823	USP <800> 60 min 30.3 min 6.2mil 		STERILE						
<b>Sterile pairs (Nitrile)</b> N8830-N8833	USP <800> 60 min 30.3 min 6.2mil 		STERILE						
<b>Extended Cuff (Nitrile)</b> N8850XPB-N8854XPB	USP <800> 6.9mil 								
<b>CHEMOPLUS™ Sterile (Nitrile)</b> CTS731S - CTS731XL	USP <800> 8.0mil 		STERILE						
<b>CHEMOPLUS™ (Nitrile)</b> CT5071G-CT5074G	USP <800> 8.0mil 								
<b>CHEMOPLUS™ (Neoprene)</b> CT21911-CT21941	USP <800> 30.4 min 9.0mil 								
<b>CHEMOPLUS™ (Latex)</b> CT0191-1-CT0194-1	USP <800> 18.0mil 								
<b>XP Plus (Latex)</b> L88HRS-L88HRXL	USP <800> 30.6 min 14.1mil 								
<b>Decontamination (Nitrile)</b> 88NDS-88NDXL	USP <800> 67.8 min 11.0mil 								

### Key

Protection against blood and bodily fluid	minimal  increased  maximum 	Cleaning 
Chemo		Smooth 
USP <800> compliant	USP <800> <sup>†</sup>	Fingertip texture 
Carmustine minimum breakthrough detection time in minutes		Full finger texture 
Thiotepa minimum breakthrough detection time in minutes		Palm and finger texture 
		Fully textured 

§ This icon indicates that the product has been recommended for chemotherapy administration given the thickness and cuff length. They also meet the ASTM D6978 standard. Please note, the majority of our products also meet this standard.

† Meets the ASTM D6978 standard  
This Cardinal Health product helps healthcare facilities comply with the USP <800> standard for the safe handling of hazardous drugs.



[cardinalhealth.com/examgloves](https://cardinalhealth.com/examgloves)

© 2019 Cardinal Health. All Rights Reserved.  
CARDINAL HEALTH, the Cardinal Health LOGO, FLEXAL, CHEMOPLUS, ESTEEM and ESSENTIAL TO CARE, are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Lit. No. 2MP19-936988 (04/2019)

This document was developed by a group of clinicians. It reflects current best practices of exam glove usage per application.